


# Vera180™

Synthetic Absorbable Lacrimal Plugs

**Vera180™** synthetic absorbable lacrimal plugs provide approximately 6 months of treatment for dry eye syndrome.

- Provides extended temporary occlusion lasting approximately **180 days**
- Ideal for treatment of **post-surgical dry eye**, and to determine potential effectiveness of a non-absorbable punctal plug
- Also indicated for use in treating the dry eye components of various **ocular surface diseases**, including contact lens intolerance



**Months**  
OF RELIEF FROM DRY EYE

## Vera180™ SINGLE PAIR

Synthetic absorbable lacrimal plugs made of poly-p-dioxanone (PDO). Absorbs in approximately 180 days. Two plugs per box.

### SIZE

0.2mm	ESK-1002S
0.3mm	ESK-1003S
0.4mm	ESK-1004S
0.5mm	ESK-1005S

## Vera180™ TEN PAIRS

Synthetic absorbable lacrimal plugs made of poly-p-dioxanone (PDO). Absorbs in approximately 180 days. 20 plugs per box (10 packs of 2 plugs each).

### SIZE

0.2mm	ESK-1002
0.3mm	ESK-1003
0.4mm	ESK-1004
0.5mm	ESK-1005