

PRODUCT COMPARISON



premiumLine

UFSK
INTERNATIONAL





surgiForce
Item No.: 2600 0002

surgiForce light
Item No.: 2600 0005

Status:
10/2021



chassis	aluminium cast technology - maximum stability and longevity	
chassis / dimensions	660 x 715 mm (incl. brake lever and foot switch)	660 x 710 mm (incl. brake lever and foot switch)
castor	4x 100 mm Ø 1x 65 mm Ø	
brake system	central brake - 99% braking power on 4 castors	
battery	strong & environmentally accumulators	
mains operation	Yes	
footswitch	 height adjustment arms & seat	 height adjustment seat
top frame	stainless steel	zinc undercoated, powder-coated
seating systems	Item No.: 2600 0004 Item No.: 2600 0002	Item No.: 2600 0006 Item No.: 2600 0005
height adjustment	Comfort Sitz 600 - 800 mm Dynamic Sitz 620 - 820 mm	Comfort Sitz 600 - 800 mm Dynamic Sitz 620 - 820 mm
seat slide	Yes (standard)	Yes (optional) Item No.: 3351 1013
armrests	stainless steel, PUR foam, sterile one-hand release 3-D adjustments - maximum flexibility	
height adjustment armrests	electrically, synchronous	gas spring provided, asynchronous
accessory armrest pair for neurosurgery	Yes (optional) Item No.: 3351 1010	No
armrest pair with linear guide	Yes (optional) Item No.: 3351 1008	
footswitch plate adapter	Yes (optional) Item No.: 3317 0000	
max. load	150 kg	

