



PRODUCT COMPARISON

performanceLine / functionalClass

UFSK
INTERNATIONAL



ophthalmic surgery





IVOM treatments



implantology maxillofacial surgery



Status:
01/2020

		500 XLE Item No.: 1500 0006	500 XLS Item No.: 1500 0005	
	chassis (LxW)	785 x 720 mm	1020 x 710 mm	
	material	zinc undercoated, powder-coated welded steel frame	zinc undercoated, powder-coated welded steel frame	
	castors	4 big double castors 4x 125 mm Ø, 360°, braking, neutral, steering, electrically conducting	4 castors 2x double castors 100 mm Ø, 360°, braking, 2x 150 mm Ø, fixed steering	
	braking system	manual		
	central brake	No, brake left/right	No, brake front castors	
	brake lever	left/right foot lever		
	steering	Yes (manual)	Yes	
	foot switch	removable foot control		
	foot switch functions	on/off (POWER) adjusting head, height	on/off (POWER) adjusting back, seat, height 3-memo- / reset-position	
	top frame	zinc undercoated, powder-coated welded steel frame	zinc undercoated, painted steel frame	
	bed surface (width)	610 mm	600 mm	
	norm rails	Yes (optional) Item No.: 3351 0047	Yes (optional) Item No.: 3321 0000 / 3320 0000	
	upholstery	back 72 mm seat 75 mm leg 60 mm	„Tundra skai“ back 100 mm seat 100 mm leg 80 mm	
	headrest	foam PUR, seamless, with hollow and neck support 	U shape, seams, small hollow 	
	headpart adjustments	electrical, manual	manual	
	pushing handles	foot end, removable, ergonomic shape stainless steel with rubber	foot end, fixed, stainless steel	



PRODUCT COMPARISON

performanceLine / functionalClass

UFSK
INTERNATIONAL

OSYS



ophthalmic surgery



IVOM treatments



implantology maxillofacial surgery



500 XLE

Item No.: 1500 0006

500 XLS

Item No.: 1500 0005

Status:
01/2020

