



PRODUCT COMPARISON

premiumLine / performanceLine

UFSK
INTERNATIONAL




ophthalmic surgery



IVOM treatments



implantology maxillofacial surgery



	eyeForce Item No.: 1800 0000	500 XLE Item No.: 1500 0006	Status: 01/2020
chassis (LxW)	900 x 662 mm	785 x 720 mm	
material	light weight chassis aluminium cast technology	zinc undercoated, powder-coated welded steel frame	
castors	4 big double castors 2x 150 mm Ø / 2x 100 mm Ø, 2x braking, neutral, steering, electrically conducting	4 big double castors 4x 125 mm Ø, 360°, braking, neutral, steering, electrically conducting	
braking system	electrically	manual	
central brake	Yes	No, brake left/right	
brake lever	electrically	left/right foot lever	
steering	Yes (electrically)	Yes (manual)	
foot switch	retractable joystick	removable foot control	
foot switch functions	on/off (POWER) adjusting head, back, seat, height	on/off (POWER) adjusting head, height	
top frame	stainless steel frame lasercut	zinc undercoated, powder-coated welded steel frame	
SES – surface extension system	Yes	No	
bed surface (width)	650 mm	610 mm	
norm rails	Yes (standard)	Yes (optional) Item No.: 3351 0047	
upholstery	back 60 mm seat 60 mm leg 50 mm	„Tundra skai“ back 72 mm seat 75 mm leg 60 mm	
headrest	foam PUR, seamless, with hollow and neck support		
headpart adjustments	electrically, manual, 2 levers (lever open = automatic adjustment)		
pushing handles	head and foot end, tiltable, ergonomic shape, joysticks included	foot end, removable, ergonomic shape, stainless steel with rubber	



PRODUCT COMPARISON

premiumLine / performanceLine

UFSK
INTERNATIONAL



ophthalmic surgery



IVOM treatments



implantology maxillofacial surgery




eyeForce

Item No.: 1800 0000

500 XLE

Item No.: 1500 0006

Status:
01/2020



Hand control functions	on/off (POWER) adjusting head, back, seat, foot, height	
battery display LED	Yes	
brake activ LED	Yes	No
steering guide activ LED	Yes	No
reset position	Yes	
memo positions	5	8
CPR button	No	Yes
max. weight	250 kg	
battery	Item No.: 3339 0000	back up battery Item No.: 3047 0000
mains operation	No	
surgery cycle	aprox. 100 without e-drive aprox. 80 with e-drive	aprox. 100
e-drive	Yes (optional) Item No.: 3316 0000	No
entry height	580 mm	standard 580 mm / 500 mm (incl. legrest extension Item No.: 3351 1028) comfort 610 mm / 530 mm (incl. legrest extension Item No.: 3351 1028)