



surgiLine

SURGEON'S CHAIR



ENJOY QUALITY AND DESIGN

surgiLine




Item No.: 2200 0008 with Formed seat (Illustration)
Item No.: 2200 0013 with Comfort seat

The quintessential timeless surgeon chair for microsurgery. The *surgiLine* is a tried and tested classic and a reliable partner in the operating room. It has a rugged steel frame yet is easy to maneuver thanks to its high-grade double castors. The seat and backrest are height-adjustable and promote a healthy sitting posture thanks to *surgiLine's* special lumbar support. The armrests can be individually adjusted for each surgeon. Available in two seating system types: Formed seat or Comfort.



Seating systems



Formed seat

Seat angle adjustment



Back angle adjustment



Back height adjustment



Upholstery / colour card

Standard colour
anthracite
Item No.: 3037 0052



Formed seat

"Tundra skai" Page 108


"Pandoria Plus" Page 109

Standard colour
black
Item No.: 3037 0235



Comfort

"Periline" Page 111



Comfort


Seat angle adjustment



Back angle adjustment



Back height adjustment



Armrests

- Stable, ergonomically shaped armrests with ample support surface (350 x 120 mm)
- Height-adjustable and equipped with a ball joint for 360-degree rotation
- Designed for personalized adjustment to the surgeon
- Optionally available with linear guide and one-hand release



Chassis

- Powder-coated tubular steel frame with sanitary, scratch-resistant high-sheen plastic cover panels
- Central braking system effective on 4 castors
- Integrated foot keypad panel for seat height adjustment
- Power supply-independent battery box with emergency stop button



Technical specifications

Lifting column

- Drive 3,000 N (IPX4)
- Lifting height 200 mm

Chassis

- L x W (incl. braking system/foot keypad) 640 x 600 mm
- Castors ø 100 mm
- Central brake acting on all four castors
- Emergency stop button

Battery unit

- Battery box with 2 batteries (lead-gel), pollution free, maintenance-free, cadmium-free 2 x 12 V, 7.2 Ah
- Universal battery charger: Input: 100 – 240 V
Output: 24 V

Seating system: Formed seat (Comfort)

- Min. sitting height 530 mm (550 mm)
- Max. sitting height 730 mm (750 mm)
- Seat angle adjustments +6°/-5° (+5°/-5°)
- Back angle adjustments +5°/-18° (+5°/-25°)

Weight 45 kg

Max. load 150 kg

Subject to technical modifications.

Technical dimensions

Seating system Formed seat



Seating system Comfort



Components



surgiTrend
Dynamic seat Item No.: 2200 0010



surgiTrend
Comfort seat Item No.: 2200 0012



surgiLine
Formed seat Item No.: 2200 0008



surgiLine
Comfort seat Item No.: 2200 0013



Backup battery

Item No.: 3047 0000



surgiTrend
Footswitch plate
480 x 715 mm, height 53 mm

Item No.: 3351 1011



surgiLine
Footswitch plate
480 x 715 mm, height 53 mm

Item No.: 3351 0014



Armrest pair with
linear guide and
one-hand release,
factory-fitted

Item No.: 3351 006



Armrest pair with
one-hand release,
factory-fitted

Item No.: 3351 1004

"Tundra skai"








High quality skai® material with fine calf leather grain and subdued matte-gloss-effect. The very extensive and trend-orientated colour spectrum enables many applications and combinations.

Product Specification

- Long-lasting
- Tear resistant
- Easy to clean
- Abrasion resistant (DIN EN ISO 6570-2)
- Lightfast
- Disinfectant resistant
- Fulfills medical standard (DIN EN ISO 10993-5+10)
- Blood- and urine resistant
- Phthalate-free
- Vegan

eyeForce	1800 0000
500 XLE	1500 0006
500 XLE comfort	1500 0006
	+ 3037 0080
500 bariatric	1500 0002
500 XLS	1500 0005
500 ECO	1500 0008
400 XL	1400 0001
surgiForce	2600 0008
surgiForce light	2600 0009
surgiLine	2200 0008
assistLine	2400 0015
motionLine	1700 0000
motionLine-E	1700 0001

+ Item No.: 3037 0...

 F6461198 Item No.:134	 F6461210 Item No.:057	 F6461636 Item No.:094	 F6461230 Item No.:058	 F6461610 Item No.:153	 F6461495 Item No.:067	 F6461456 Item No.:001	 F6461338 Item No.:008
 F6461288 Item No.:007	 F6461550 Item No.:070	 F6461608 Item No.:075	 F6461496 Item No.:068	 F6461609 Item No.:076	 F6461581 Item No.:073	 F6461196 Item No.:051	 F6461453 Item No.:062
 F6461356 Item No.:059	 F6461381 Item No.:119	 F6461354 Item No.:011	 F6461350 Item No.:010	 F6461199 Item No.:054	 F6461454 Item No.:063	 F6461557 Item No.:072	 F6461556 Item No.:071
 F6461458 Item No.:066	 F6461457 Item No.:065	 F6461607 Item No.:074	 F6461542 Item No.:069	 F6461343 Item No.:009	 F6461543 Item No.:109	 F6461455 Item No.:064	 F6461273 Item No.:006
 F6461205 Item No.:056	 F6461204 Item No.:055	 F6461357 Item No.:060	 F6461195 Item No.:050	 F6461197 Item No.:052			

The colours displayed on your device may differ slightly from the original colours. If you wish to receive original colour charts, please contact us at info@ufsk-osys.com.

"Pandoria Plus"



skai® Pandoria Plus is a high-quality upholstery material, whose fine micropore structure gives a certain unobtrusiveness, which is refined by an attractive tone-in-tone effect. The very pleasant feel of the material and the excellent material characteristic and the broad, modern color spectrum make skai® Pandoria Plus universally usable.


Product Specification

- Long-lasting
- Tear resistant
- Easy to clean
- Abrasion resistant (DIN EN ISO 5470-2)
- Resistant to oil and fat
- Lightfast
- Disinfectant resistant
- Flame protection
- Fulfills medical standard
- Resistant to sudor, blood and urine
- 0161-21 (IMO FTPC Part 8)
- Antibacterial according to JIS Z 2801
- Vegan


<i>eyeForce</i>		1800 0000
<i>500 XLE</i>		1500 0006
<i>500 XLE comfort</i>		1500 0006
		+ 3037 0080
<i>500 bariatric</i>		1500 0002
<i>500 XLS</i>		1500 0005
<i>500 ECO</i>		1500 0008
<i>400 XL</i>		1400 0001
<i>surgiForce</i>	Formed seat	2600 0008
<i>surgiForce light</i>	Formed seat	2600 0009
<i>surgiLine</i>	Formed seat	2200 0008
<i>assistLine</i>		2400 0015
<i>motionLine</i>		1700 0000
<i>motionLine-E</i>		1700 0001

+ Item No.: 3037 0136

+ Item No.: 3037 0...




weiß F6413046
Item No.:053



hellbeige F6413047
Item No.:141




kiesel F6413048
Item No.:158



birke F6413049
Item No.:159



smoke F6413050
Item No.:221



silbergrau F6413075
Item No.:275



platin F6413053
Item No.:161



anthrazit F6413076
Item No.:276



schwarz F6413052
Item No.:170



schoko F6413051
Item No.:160



gletscher F6413055
Item No.:162



sky F6413074
Item No.:277



royal F6413056
Item No.:163



smaragd F6413073
Item No.:278



kaktus F6413077
Item No.:281



feuer F6413058
Item No.:166



kirsche F6413079
Item No.:279



fuchsia F6413078
Item No.:280



mais F6413057
Item No.:164



limone F6413054
Item No.:165

The colours displayed on your device may differ slightly from the original colours. If you wish to receive original colour charts, please contact us at info@ufsk-osys.com.

"Periline"



Product Specification

- Without Phthalates
- Biocompatibility:
Cytotoxicity (EN ISO 10 993-5) and
test for irritation and delayed-type
hypersensitivity (EN ISO 10 993-10)
- Urine resistance
- Blood resistance
- Resistance to perspiration
- Resistance to spitting
- Light fastness
- Color fastness to rubbing
- Flame retardant

<i>surgiForce</i>	Comfort seat	2600 0004
<i>surgiForce light</i>	Comfort seat	2600 0006
<i>surgiTrend</i>	Comfort seat	2200 0012
<i>surgiLine</i>	Comfort seat	2200 0013
<i>assistPro F</i>		2400 0005
<i>assistPro H</i>		2400 0009
<i>easyPro F</i>		2400 0010
<i>easyPro H</i>		2400 0006
<i>easyPro R</i>		2400 0013
<i>easyLine F</i>		2400 0020
<i>easyLine H</i>		2400 0021
<i>easyLine R</i>		2400 0022
<i>easyLine Hs</i>		2400 0023
<i>easyLine Rs</i>		2400 0024

+ Item No.: 3037 0...

				
kastanie 015817 Item No.:240	light brown 016037 Item No.:239	creme 015812 Item No.:225	weiß 010881 Item No.:224	stein 015813 Item No.:236
				
saphir 015814 Item No.:231	delft 010766 Item No.:241	antistatisch schwarz Item No.:223	schwarz 011499 Item No.:235	anthracit 016046 Item No.:234
				
azur 015815 Item No.:230	aubergine 016041 Item No.:229	lavendel 016042 Item No.:228	rot 010878 Item No.:238	orange 016040 Item No.:237
				
kiwi 015816 Item No.:233	lime green 016043 Item No.:232	lemon 016039 Item No.:226	safran 010879 Item No.:227	

The colours displayed on your device may differ slightly from the original colours. If you wish to receive original colour charts, please contact us at info@ufsk-osys.com.

Upgrade your OR to perfection.



eyeForce



500 XLE



500 XLE comfort



500 bariatric



500 XLS



500 ECO



400 XL



motionLine



motionLine-E



surgiForce



surgiForce light



surgiTrend



surgiLine



stayClean wipes



stayClean



staySoft



stayClean roll



supportLine



supportLine plus



assistTrend



assistPro F



assistPro H



easyPro F



easyPro H



easyPro R



easyLine F



easyLine H



easyLine R



easyLine Hs



easyLine Rs



assistLine

CONTACT US



UFSK
INTERNATIONAL

OSYS

UFSK-International OSYS GmbH
Kirchhoffstraße 1
93055 Regensburg
Germany

Phone +49 (0) 941 78862 15
Fax +49 (0) 941 78862 35
Mail info@ufsk-osys.com
Web ufsk-osys.com

MADE IN GERMANY
Design | Construction | Manufacturing