

Linear guides

To give you that extra margin of flexibility:
Armrests with linear guides! Applying pressure with your forearm will make the arms "glide" effortlessly in longitudinal direction while letting go will stop them by self-locking mechanism - available as an option.



Armrests

Made of molded PUR foam - pleasant to the touch, no slick leatherette cover - for guaranteed secure hold. Suitable for disinfectants.



Upholstery / colour card

Standard colour
black
Item No.: 3037 0235



"Periline" Page 111

Comfort

Standard colour
anthracite
Item No.: 3037 0052



"Tundra skai" Page 108
"Pandoria Plus" Page 109

Formed seat

Technical specifications

Lifting column

- Lifting height of the seat 200 mm
- Lifting height of the armrests 130 mm
- Electromotive adjustment via the foot control elements

Dimensions

- Max. LxWxH (incl. armrests) 860x750x1,270 mm

Chassis

- Castors ø 100 mm
- 5th castor ø 65 mm
- Central brake acting on all four castors
- Emergency stop button

Battery unit

- Lead-gel battery, 2 pcs. 12 V, 2.9 A
- Universal battery charger: Input: 100 – 240 V
Output: 24 V

Sitzsysteme: Comfort (Formed seat)

- Min. sitting height 600 mm (580 mm)
- Max. sitting height 800 mm (780 mm)
- Seat position adjustment (forward/back) 120 mm
- Seat angle adjustment -6° (-14°)
- Back angle adjustment -10° (-19° / +5°)
- Seat adjustment mechanism 113 mm (113 mm)

Weight

74 kg


Max. load

150 kg

Techn. Änderungen vorbehalten.

Technical dimensions

Seating system Comfort



Seating system Formed seat

