



supportLine

INSTRUMENT TABLE



ENJOY QUALITY AND DESIGN

E

supportLine

Item No.: 3351 0026

The *supportLine* mobile instrument table is designed to be an ergonomic aid to supplement treatment and surgeon's chairs.

**The practical
little helpers**



Technical specifications

Chassis

- Footprint 680x400 mm
- 4 castors ø 50 mm

Lifting column

- Height adjustment 730 – 985 mm

Stainless steel tray

660x540 mm

Upper frame

715x540 mm

Weight

17,5 kg

Max. load

10 kg

Subject to technical modifications.

The perfect addition for your OR

Pneumatic spring technology/foot pedal

The integrated pneumatic spring technology enables the instrument tray to be stopped and positioned smoothly at a suitable working height for the surgeon. Infinitely variable height adjustment via foot pedal.



supportLine plus



Item No.: 3351 0045

The *supportLine*'s maneuverability, optimised stability and resistance to tipping, its stainless steel design for the user's operating area and its hard-wearing, hygienic surfaces are designed specifically for versatile use in the operating room. The *supportLine plus* mobile instrument table is a 20 cm higher version.



Technical specifications

Chassis

- Footprint 680x400 mm
- 4 castors ø 50 mm

Lifting column

- Height adjustment 930 – 1,185 mm

Stainless steel tray

660x540 mm

Upper frame

715x540 mm

Weight

17,5 kg

Max. load

10 kg

Subject to technical modifications..

Tray

The removable, high-quality stainless steel tray is heat-resistant, easy to clean and disinfect and offers plenty of space for holding surgical instruments with an area of 66 x 54 cm.



Push handle/double castors

Its wide push handle and maneuverable double castors make the *supportLine* easy to position.



Upgrade your OR to perfection.



eyeForce



500 XLE



500 XLE comfort



500 XLS



500 ECO



400 XL



motionLine



motionLine-E



surgiForce



surgiForce light



surgiTrend



surgiLine



stayClean wipes



stayClean



staySoft



stayClean roll



supportLine



supportLine plus



assistTrend



assistPro F



assistPro H



easyPro F



easyPro H



easyPro R



easyLine F



easyLine R



easyLine Rs



assistLine

CONTACT US



UFSK
INTERNATIONAL

OSYS

UFSK-International OSYS GmbH
Kirchhoffstraße 1
93055 Regensburg
Germany

Phone +49 (0) 941 78862 15
Fax +49 (0) 941 78862 35
Mail info@ufsk-osys.com
Web ufsk-osys.com

MADE IN GERMANY
Design | Construction | Manufacturing