

# Refraction unit ELITE



ELITE



## Refraction unit ELITE

- Motorized sliding table top for two instruments;
- Elevating table top (25 cm);
- Elevating patient chair (15 cm);
- Removable patient chair for wheelchair access;
- Column with overhead LED lamp with adjustable illumination;
- Column with support for chart projector ;
- Integrated adjustable power transformers for slit lamp/ophthalmometer;
- Chin rest support plates;
- 230V outlet for additional instrument e. g. auto lensmeter, auto phoropter, etc.;
- Head rest for the chair;
- Cable drag chain for data cables;
- Safety switches for horizontal and vertical movement of the table.


## Optional

- Motorized phoropter arm
- Foot rest (either fixed or foldable)
- Foldable arm rests
- Chair with reclining back
- Manual (13cm) or motorized (14cm) chair sliding forward – backwards
- Foot switch control for up/down movement of the chair
- Installed data cables (Fire Wire, USB, RS232, etc.)
- Pillow for children
- Doctor's desk according customer specifications
- Left-handed version

Refraction unit

**Med**  
instrus

# Refraction unit ELITE




1. Cabinet on castors with 4 drawers. 2. Column with overhead halogen lamp with adjustable illumination. 3. Chair control. 4. Control panel 5. Elevating table top. 6. Additional shelf.



## Technical information

Voltage	~230V 50/60Hz
Total weight	210 kg
Max height	186 cm
Table height from the ground	From 73cm to 98cm (25cm lift)
Chair height from the ground*	From 48cm to 63cm (15cm lift)
Chair rotation*	360°
Overhead LED lamp	2x12V 5W
Instruments power supply	~230V, ~12V adj., ~6V adj.

\*Can vary depending on options.



UAB „Medinstrus“  
Respublikos 20, LT - 89224 Mažeikiai,  
Lietuva (Lithuania).

Tel: +370 443 25451  
Fax: +370 443 25455  
E-mail: info@medinstrus.lt  
Web: www.medinstrus.lt

**Med**  
**instrus**